

ZeniABSORB

Super Absorbent Dressing

ZeniMedical



Super Absorbent Dressing

ZeniABSORB is a sterile dressing for use on heavily exuding wounds, consisting of an oxygen permeable, waterproof nonwoven layer, a super absorbent pad layer, hydro-penetrating nonwoven layer and silicone nonadherent contact layer. The layered construction facilitates fluid management to provide optimal moist wound environment and locks in exudate.

Product Features

- Super Absorbency
- Non-Adherent to Wound Bed
- Breathable, Waterproof Protective Backing Prevents Strike Through
- Available in Multiple Sizes
- Protects Against Maceration
- Can be Repositioned upon Application
- Gentle on Sensitive and Fragile Skin
- Used as Primary or Secondary Dressing
- Highly Conformable and Comfortable
- Easy Application

Indications for Use

- Pressure Ulcers
- Diabetic Ulcers
- Moderate to Heavily Draining Wounds
- Leg Ulcers
- Surgical Wounds

Contraindications

ZeniABSORB is contraindicated for use on third degree burns, dry or low exuding wounds, eyes, mucous membranes, as an internal dressing/sponge or wound cavities due to the dressing swells during absorption.

Item Number	Description	Size	Packaging	HCPCS Code
80045	ZeniABSORB	4" x 5"	10/Box	Pending
80069	ZeniABSORB	6" x 9"	10/Box	Pending

* Providers are responsible for determining appropriate billing codes for Medicare Program Claims. This information is included for convenience and is not intended, no should be considered billing advice.

To order, call: 888.945.1492 • Fax: 888.945.1492 • E-mail: orders@Zenimedical.com • Knoxville, TN 37909 • ZeniMedical.com